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1. Bergstrom A, Skov T.H., Bahl M.L., Roager H.M., Christensen L.B., Ejlerskov K.T. et al. Establishment of intestinal microbiota during early life: a longitudinal explorative study of a large cohort of Danish infants. *Appl. Environ. Microbiol.* 2014; 80(9): 2889–900.
2. Arrieta M, et al. Early infancy microbial and metabolic alterations affect risk of childhood asthma. *Sci Transl Med.* 2015; 7:1–14.
3. Fujimura K.E., Sitarik A.R., Havstad S, Lin D.L., Levan S, Fadrosch D, Panzer A.R., LaMere B, Elze Rackaityte, Nicholas W Lukacs, Ganesa Wegienka, Homer A Boushey, Dennis R Ownby, Edward M Zoratti, Albert M Levin, Christine C Johnson, Susan V Lynch. Neonatal gut microbiota associates with childhood multisensitized atopy and T cell differentiation. *Nat Med.* 2016 Oct;22(10):1187-1191. doi: 10.1038/nm.4176.
4. Sidhu M, van der Poorten D. The gut microbiome. *Aust Fam Physician* 2017.
5. Bäckhed F, Fraser Claire M, Ringel Y, Sanders M.E., Sartor R.B., Sherman P.M., et al. Defining a healthy human gut microbiome: current concepts, future directions, and clinical applications. *Cell Host Microbe* 2012.
6. Hollister E.B., Gao C, Versalovic J. Compositional and functional features of the gastrointestinal microbiome and their effects on human health. *Gastroenterology.* 2014.
7. Guarner, F.; Malagelada J.R. Gut flora in health and disease. *Lancet* 2003, 361, 512–519.
8. Stecher, B.; Hardt W. The role of microbiota in infectious disease. *Trends Microbiol.* 2008, 16, 107-114.
9. Machiels K., Joossens M., Sabino J., de Preter V., Arijis I., Eeckhaut V., Ballet V., Claes K., van Immerseel F., Verbeke K. et al. A decrease of the butyrate-producing species *Roseburia hominis* and *Faecalibacterium prausnitzii* defines dysbiosis in patients with ulcerative colitis. *Gut* 2014, 63, 1275–1283.
10. McLroy J, Ianiro G, Mukhopadhyaya I, Hansen R, Hold GL. Review article: the gut microbiome in inflammatory bowel disease – avenues for microbial management. *Aliment Pharmacol Ther.* 2018.
11. Zhang Y-J, Li S, Gan R-Y, Zhou T, Xu D-P, Li H-B. Impacts of gut bacteria on human health and diseases. *Int J Mol Sci.* 2015.
12. Tang M.S., Poles J, Leung J.M., Wolff M.J., Davenport M, Lee S.C., Lim Y.A., Chua K.H., Loke P, Cho I. Inferred metagenomic comparison of mucosal and fecal microbiota from individuals undergoing routine screening colonoscopy reveals similar differences observed during active inflammation. *Gut Microbes.* 2015;6(1):48-56. doi: 10.1080/19490976.2014.1000080.
13. McCoy A.N. et al. *Fusobacterium* is associated with colorectal adenomas. *PLOS ONE* 8, e53653 (2013).
14. Loftus M, Hassouneh SA-D, Yooseph S. Bacterial community structure alterations within the colorectal cancer gut microbiome *BMC Microbiol.* 2021 Mar 31;21(1):98. doi: 10.1186/s12866-021-02153-x.
15. Osman M.A., Neoh H-M, Ab Mutalib N-S, Chin S.F., Mazlan L, Ali RAR, Zakaria A.D., Ngiu C.S., Ang M.Y., Jamal R. *Parvimonas micra*, *Peptostreptococcus stomatis*, *Fusobacterium nucleatum* and *Akkermansia muciniphila* as a four-bacteria biomarker panel of colorectal cancer. *Sci Rep.* 2021 Feb 3;11(1):2925. doi: 10.1038/s41598-021-82465-0.
16. Moran, C.P.; Shanahan, F. Gut microbiota and obesity: Role in aetiology and potential therapeutic target. *Best Pract. Res. Clin. Gastroenterol.* 2014, 28, 585–597.
17. Crovesy, L.; Masterson, D.; Rosado, E.L. Profile of the gut microbiota of adults with obesity: A systematic review. *Eur. J. Clin.Nutr.* 2020.
18. Atzeni A, Galié S, Muralidharan J, Babio N, Tinahones FJ, Vioque J, Corella D, Castañer O, Vidal J, Moreno-Indias I, Torres-Collado L, Fernández-Carrión R, Fitó M, Olbeyra R, Martínez-González MA, Bulló M, Salas-Salvadó. Gut microbiota profile and changes in body weight in elderly subjects with overweight/obesity and metabolic syndrome. *J. Microorganisms.* 2021 Feb 10;9(2):346. doi: 10.3390/microorganisms9020346.
19. Chen Y, Zhou J, Wang L. Role and Mechanism of Gut Microbiota in Human Disease. *Front Cell Infect Microbiol.* 2021 Mar 17;11:625913. doi: 10.3389/fcimb.2021.625913.

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4. The Role of FODMAPs in Irritable Bowel Syndrome. Susan J Shepherd 1, Emma Halmos, Simon Glance.

5. Salminen S, Isolauri E. Intestinal colonization, microbiota and probiotics. *J Pediatr.* 2006;149:S115–S120.

6. Hébuterne X. Gut changes attributed to ageing: effects on intestinal microflora. *Curr Opin Clin Nutr Metab Care.* 2003;6:49–54.

7. Howard M, Gordon D, Garleb K, Kerley M. Dietary fructooligosaccharide, xylooligosaccharide and gum arabic have variable effects on cecal and colonic microbiota and epithelial cell proliferation in mice and rats. *J Nutr.* 1995;

8. Hernández-Hernández O, Marín-Manzano M, Rubio L, Moreno F, Sanz M, Clemente A. Monomer and linkage type of galacto-oligosaccharides affect their resistance to ileal digestion and prebiotic properties in rats. *J Nutr.* 2012.

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
1. <http://www.bacterio.net/lactobacillus.html>

2. Dheeney D., Gareau M.G., Marco M.L. Intestinal *Lactobacillus* in health and disease, a driver or just along for the ride? *Curr. Opin. Biotechnol.* 2018;49:140–147. doi: 10.1016/j.copbio.2017.08.004.

3. Amadoro C., Rossi F., Pallotta M.L., Gasperi M., Colavita G. Traditional dairy products can supply beneficial microorganisms able to survive in the gastrointestinal tract. *LWT Food Sci. Technol.* 2018.

4. Walter J. Ecological role of lactobacilli in the gastrointestinal tract: Implications for fundamental and biomedical research. *Appl. Environ. Microbiol.* 2008.

5. Liu Y.W., Su Y.W., Ong W.K., Cheng T.H., Tsai Y.C. Oral administration of *Lactobacillus plantarum* K68 ameliorates DSS-induced ulcerative colitis in BALB/c mice via the anti-inflammatory and immunomodulatory activities. *Int. Immunopharmacol.* 2011.

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1. Bacteroides Thetaiotaomicron. Nathan T Porter , Ana S Luis , Eric C Martens.

2. Hooper L.V., Wong M.H., Thelin A, et al. Molecular analysis of commensal host-microbial relationships in the intestine. Science. 2001.

3. Wrzosek L, Miquel S, Noordine M.L., et al. . Bacteroides thetaiotaomicron and faecalibacterium prausnitzii influence the production of mucus glycans and the development of goblet cells in the colonic epithelium of a gnotobiotic model rodent. BMC Biol. 2013.

4. Kelly D, Campbell J.I., King T.P., et al. . Commensal anaerobic gut bacteria attenuate inflammation by regulating nuclear-cytoplasmic shuttling of PPAR-gamma and rela. Nat Immunol. 2004.

5. Bacteroides thetaiotaomicron Ameliorates Colon Inflammation in Preclinical Models of Crohn's Disease. Delday M, Mulder I, Logan E.T., Grant G. Inflamm Bowel Dis. 2019.

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3. De Vos W.M. Microbe Profile: *Akkermansia muciniphila*: A conserved intestinal symbiont that acts as the gatekeeper of our mucosa. *Microbiology.* 2017.

4. Zhang X., Shen D., Fang Z., Jie Z., Qiu X., Zhang C., Chen Y., Ji L. Human gut microbiota changes reveal the progression of glucose intolerance. *PLoS ONE.* 2013.

5. Swidsinski A., Dorffel Y., Loening-Baucke V., Theissig F., Ruckert J.C., Ismail M., Rau W.A., Gaschler D., Weizenegger M., Kuhn S., et al. Acute appendicitis is characterised by local invasion with *Fusobacterium nucleatum/necrophorum*. *Gut.* 2011.

6. Karlsson C.L., Onnerfalt J., Xu J., Molin G., Ahre S., Thorngren-Jerneck K. The microbiota of the gut in preschool children with normal and excessive body weight. *Obesity.* 2012.

7. Dao M.C., Everard A., Aron-Wisniewsky J., Sokolovska N., Prifti E., Verger E.O., Kayser B.D., Levenez F., Chilloux J., Hoyles L., et al. *Akkermansia muciniphila* and improved metabolic health during a dietary intervention in obesity: Relationship with gut microbiome richness and ecology. *Gut.* 2015.

8. Zhang X., Shen D., Fang Z., Jie Z., Qiu X., Zhang C., Chen Y., Ji L. Human gut microbiota changes reveal the progression of glucose intolerance. *PLoS ONE.* 2013.

9. Clarke S.F., Murphy E.F., O'Sullivan O., Lucey A.J., Humphreys M., Hogan A., Hayes P., O'Reilly M., Jeffery I.B., Wood-Martin R., et al. Exercise and associated dietary extremes impact on gut microbial diversity. *Gut.* 2014.

10. Barton W., Penney N.C., Cronin O., Garcia-Perez I., Molloy M.G., Holmes E., Shanahan F., Cotter P.D., O'Sullivan O. The microbiome of professional athletes differs from that of more sedentary subjects in composition and particularly at the functional metabolic level. *Gut.* 2018.

11. Jangi S., Gandhi R., Cox L.M., Li N., von Glehn F., Yan R., Patel B., Mazzola M.A., Liu S., Glanz B.L., Cook S., Tankou S., Stuart F., Melo K., Nejad P., Smith K., Topcuolu B.D., Holden J., Kivisäkk P., Chitnis T., De Jager P.L., Quintana F.J., Gerber G.K., Bry L., Weiner H.L. Alterations of the human gut microbiome in multiple sclerosis. *Nat Commun.* 2016, 7:12015. doi: 10.1038/ncomms12015.

12. Cekanaviciute E., Yoo B.B., Runia T.F., Debelius J.W., Singh S., Nelson C.A., Kanner R., Bencosme Y., Lee Y.K., Hauser S.L., Crabtree-Hartman E., Sand I.K., Gacias M., Zhu Y., Casaccia P., Cree B.A.C., Knight R., Mazmanian S.K., Baranzini S.E. Gut bacteria from multiple sclerosis patients modulate human T cells and exacerbate symptoms in mouse models. *Proc Natl Acad Sci U S A.* 2017 Oct 3;114(40):10713-10718. doi: 10.1073/pnas.1711235114.

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
1. Oliveira F.A., Paludo K.S., Arend L.N., Farah S.M., Pedrosa F.O., Souza E.M., et al. Virulence characteristics and antimicrobial susceptibility of uropathogenic *Escherichia coli* strains. *Genet Mol Res.* 2011;10:4114–25.
2. Kim H.J., Koo M, Jeong A.R., Baek S.Y., Cho J.I., Lee S.H., et al. Occurrence of pathogenic *Escherichia coli* in commercially available fresh vegetable products in Korea. *J Korean Soc Appl Biol Chem.* 2014;57:367–70.
3. Croxen M.A., Law R.J., Scholz R, Keeney K.M., Wlodarska M, Finlay B.B., et al. Recent advances in understanding enteric pathogenic *Escherichia coli*. *Clin Microbiol Rev.* 2013;26:822–80.
4. Kaper J.B., Nataro J.P., Mobley H.L. 2004. Pathogenic *Escherichia coli*. *Nat. Rev. Microbiol.* 2:123–140.
5. Croxen M.A., Finlay B.B. 2010. Molecular mechanisms of *Escherichia coli* pathogenicity. *Nat. Rev. Microbiol.* 8:26–38.
6. World Health Organization. 2012. World health statistics 2012. WHO Press, Geneva, Switzerland.
7. Nataro J.P., Kaper J.B. 1998. Diarrheagenic *Escherichia coli*. *Clin. Microbiol. Rev.* 11:142–201.
8. Nataro J.P., Kaper J.B. 1998. Diarrheagenic *Escherichia coli*. *Clin. Microbiol. Rev.* 11:142–201.
9. Wang L, Rothmund D, Curd H, Reeves P.R. 2003. Species-wide variation in the *Escherichia coli* flagellin (H-antigen) gene. *J. Bacteriol.* 185:2936–2943.
10. DebRoy C, Roberts E, Fratamico P.M. 2011. Detection of O antigens in *Escherichia coli*. *Anim. Health Res. Rev.* 12:169–185
11. Swaminathan B, Barrett T.J., Hunter S.B., Tauxe RV. 2001. PulseNet: the molecular subtyping network for foodborne bacterial disease surveillance, United States. *Emerg. Infect. Dis.* 7:382–389.
12. Существующие открытые базы данных по MLST – <https://cge.cbs.dtu.dk/services/MLST>/<https://bigsd.pasteur.fr/https://pubmlst.org/>
13. Donnenberg M.S., Tacket C.O., James S.P., Losonsky G, Nataro J.P., Wasserman S.S., Kaper J.B., Levine M.M. 1993. Role of the *eaeA* gene in experimental enteropathogenic *Escherichia coli* infection. *J. Clin. Invest.* 92:1412–1417.
14. Nataro J.P., Kaper J.B. 1998. Diarrheagenic *Escherichia coli*. *Clin. Microbiol. Rev.* 11:142–201.
15. Goosney D.L., Gruenheid S., Finlay B.B. 2000. Gut feelings: enteropathogenic *E. coli*(EPEC) interactions with the host. *Annu. Rev. Cell Dev. Biol.* 16:173–189.
16. Nguyen R.N., Taylor L.S., Tauschek M, Robins-Browne R.M. 2006. Atypical enteropathogenic *Escherichia coli* infection and prolonged diarrhea in children. *Emerg. Infect. Dis.* 12:597–603.
17. Subramanian K, Selvakkumar C, Vinaykumar K.S., Goswami N, Meenakshisundaram S, Balakrishnan A, Lakshmi BS. 2009. Tackling multiple antibiotic resistance in enteropathogenic *Escherichia coli* (EPEC) clinical isolates: a diarylheptanoid from *Alpinia officinarum* shows promising antibacterial and immunomodulatory activity against EPEC and its lipopolysaccharide-induced inflammation. *Int. J. Antimicrob. Agents* 33:244–250
18. Le Bouguéne C, Servin A.L. 2006. Diffusely adherent *Escherichia coli* strains expressing Afa/Dr adhesins (Afa/Dr DAEC): hitherto unrecognized pathogens. *FEMS Microbiol. Lett.* 256:185–194.

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3. Janda J.M., Abbott S.L., Cheung W.K., Hanson D.F. (1994) Biochemical identification of citrobacteria in the clinical laboratory. *J Clin Microbiol* 32:1850–18542

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
3. Singh N.K., Bezdán D, Checinska Sielaff A, Wheeler K, Mason C.E., Ven-kateswaran K. 2018. Multi-drug resistant *Enterobacter bugandensis* species isolated from the International Space Station and comparative genomic analyses with human pathogenic strains. *BMC Microbiol* 18:175. <https://doi.org/10.1186/s12866-018-1325-2>.

4. Akbari M, Bakhshi B, Najar Peerayeh S. 2016. Particular distribution of *Enterobacter cloacae* strains isolated from urinary tract infection within clonal complexes. *Iran Biomed J* 20:49–55.

5. Bertrand X, Hocquet D, Boisson K, Siebor E, Plésiat P, Talon D. 2003. Molecular epidemiology of *Enterobacteriaceae* producing extended-spectrum-lactamase in a French university-affiliated hospital. *Int J Antimicrob Agents* 22:128–133. [https://doi.org/10.1016/S0924-8579\(03\)00098-0](https://doi.org/10.1016/S0924-8579(03)00098-0).

6. Chow J.W., Fine M.J., Shlaes D.M., Quinn J.P., Hooper D.C., Johnson M.P., Ramphal R, Wagener M.M., Miyashiro D.K., Yu V.L. 1991. *Enterobacter bacteremia*: clinical features and emergence of antibiotic resistance during therapy. *Ann Intern Med* 115:585–590. <https://doi.org/10.7326/0003-4819-115-8-585>.

7. Paauw A, Caspers M.P., Leverstein-van Hall M.A., Schuren F.H., Montijn R.C., Verhoef J, Fluit A.C. 2009. Identification of resistance and virulence factors in an epidemic *Enterobacter hormaechei* outbreak strain. *Microbiology* 155:1478–1488. <https://doi.org/10.1099/mic.0.024828-0>.

8. Morand P.C., Billoet A, Rottman M, Sivadon-Tardy V, Eyrolle L, Jeanne L, Tazi A, Anract P, Courpied J.P., Poyart C, Dumaine V. 2009. Specific distribution within the *Enterobacter cloacae* complex of strains isolated from infected orthopedic implants. *J Clin Microbiol* 47:2489–2495. <https://doi.org/10.1128/JCM.00290-09>.

9. Allerberger F, Koeuth T, Lass-Flörl C, Dierich M.P., Putensen C, Schmutzhard E, Mohsenipour I, Grundmann H, Hartung D, Bauernfeind A, Eberlein E, Lupski J.R. 1996. Epidemiology of infections due to multi-resistant *Enterobacter aerogenes* in a university hospital. *Eur J Clin Microbiol Infect Dis* 15:517–521. <https://doi.org/10.1007/BF01691323>

10. Wu W, Feng Y, Zong Z. 2018. *Enterobacter sichuanensis* sp. nov., recovered from human urine. *Int J Syst Evol Microbiol* 68:3922–3927. <https://doi.org/10.1099/ijsem.0.003089>.

11. Davin-Regli A, Pagès J.M. 2015. *Enterobacter aerogenes* and *Enterobacter cloacae*: versatile bacterial pathogens confronting antibiotic treatment. *Front Microbiol* 6:392. <https://doi.org/10.3389/fmicb.2015.00392>.

12. Davin-Regli A, Masi M, Bialek S, Nicolas-Chanoine M.H., Pagès J.M. 2016. Antimicrobial resistance and drug efflux pumps in *Enterobacter* and *Klebsiella*, p 281–306. In Li X-Z, Elkins C.A., Zgurskaya H.I. (ed), *Efflux-mediated drug resistance in bacteria: mechanisms, regulation and clinical implications*. Springer International Publishing, Basel, Switzerland.

13. Fernández-Baca V, Ballesteros F, Hervás J.A., Villalón P, Domínguez M.A., Benedí V.J., Albertí S. 2001. Molecular epidemiological typing of *Enterobacter cloacae* isolates from a neonatal intensive care unit: three-year prospective study. *J Hosp Infect* 49:173–182. <https://doi.org/10.1053/jhin.2001.1053>.

14. Fluit A.C., Schmitz F.J., Verhoef J. 2001. Multi-resistance to antimicrobial agents for the ten most frequently isolated bacterial pathogens. *Int J Antimicrob Agents* 18:147–160. [https://doi.org/10.1016/S0924-8579\(01\)00357-0](https://doi.org/10.1016/S0924-8579(01)00357-0).

15. Yurdakul, D., Yazgan-Karataş, A., & Şahin, F. (2015). *Enterobacter* Strains Might Promote Colon Cancer. *Current Microbiology*, 71(3), 403–411. doi:10.1007/s00284-015-0867-x

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3. Tiwana H, Natt R.S., Benitez Brito R, Shah S, Wilson C, Bridger S, et al. Correlation between the immune responses to collagens type I, III, IV and V and *Klebsiella pneumoniae* in patients with Crohn's disease and ankylosing spondylitis. *Rheumatology*. 2001;40(1):15-23.

4. Ebringer A, Rashid T, Tiwana H, Wilson C. A possible link between Crohn's disease and ankylosing spondylitis via *Klebsiella* infections. *Clinical Rheumatology*. 2007;26(3):289-97.

5. Rashid T, Ebringer A, Wilson C. The Role of *Klebsiella* in Crohn's Disease with a Potential for the Use of Antimicrobial Measures. *International Journal of Rheumatology*. 2013;2013:8.

6. Holden V.I., Breen P, Houle S, Dozois C.M., Bachman M.A. *Klebsiella pneumoniae* Siderophores Induce Inflammation, Bacterial Dissemination, and HIF-1 α Stabilization during Pneumonia. *mBio*. 2016;7(5):e01397-16

7. Lee I-A, Kim D-H. *Klebsiella pneumoniae* increases the risk of inflammation and colitis in a murine model of intestinal bowel disease. *Scandinavian Journal of Gastroenterology*. 2011;46(6):684-93.

8. Zhou W, Cao Q, Peng Y, Zhang Q-J, Castrillon D.H., DePinho R.A., et al. FoxO4 Inhibits NF- κ B and Protects Mice Against Colonic Injury and Inflammation. *Gastroenterology*. 2009; 137(4):1403-14.

9. Marcinkiewicz J, Cizek M, Bobek M, Strus M, Heczko P.B., Kurnyta M, et al. Differential inflammatory mediator response in vitro from murine macrophages to lactobacilli and pathogenic intestinal bacteria. *International Journal of Experimental Pathology*. 2007; 88(3):155-64.

10. Chung D.R., Lee H, Park M.H., Jung S, Chang H, Kim Y, Son J.S., Moon C, Kwon K.T., Ryu S.Y., Shin S.Y., Ko K.S., Kang C, Peck K.R., Song J: Fecal carriage of serotype K1 *Klebsiella pneumoniae* ST23 strains closely related to liver abscess isolates in Koreans living in Korea. *Eur J Clin Microbiol Infect Dis* 2012, 31:481-6.

11. Siu L.K., Fung C-P, Chang F-Y, Lee N, Yeh K-M, Koh T.H., et al. Molecular Typing and Virulence Analysis of Serotype K1 *Klebsiella pneumoniae* Strains Isolated from Liver Abscess Patients and Stool Samples from Noninfectious Subjects in Hong Kong, Singapore, and Taiwan. *Journal of Clinical Microbiology*. 2011; 49(11):3761-5.

12. Wang J.H., Liu Y.C., Lee S.S., Yen M.Y., Chen Y.S., Wang J.H., et al. Primary liver abscess due to *Klebsiella pneumoniae* in Taiwan. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 1998;26(6):1434-8.

13. Tsai F-C, Huang Y-T, Chang L-Y, Wang J-T. Pyogenic Liver Abscess as Endemic Disease, Taiwan. *Emerging Infectious Diseases*. 2008;14(10):1592-600.

14. Lai Y-C, Lin A-C, Chiang M-K, Dai Y-H, Hsu C-C, Lu M-C, et al. Genotoxic *Klebsiella pneumoniae* in Taiwan. *PLOS ONE*. 2014;9(5):e96292.

15. Broberg C.A., Palacios M, Miller V.L. *Klebsiella*: a long way to go towards understanding this enigmatic jet-setter. *F1000Prime Rep* 2014;6:64.

16. Beaugerie L, Metz M, Barbut F, et al. *Klebsiella oxytoca* as an agent of antibiotic-associated hemorrhagic colitis. *Clin Gastroenterol Hepatol* 2003;1:370-6.

17. Koga H, Aoyagi K, Yoshimura R, Kimura Y, Iida M, Fujishima M. Can quinolones cause hemorrhagic colitis of late onset? Report of three cases. *Dis Colon Rectum* 1999;

18. Sakurai Y, Tsuchiya H, Ikegami F, Funatomi T, Takasu S, Uchikoshi T. Acute right-sided hemorrhagic colitis associated with oral administration of ampicillin. 1979;

19. Higaki M, Chida T, Takano H, Nakaya R. Cytotoxic component(s) of *Klebsiella oxytoca* on HEP-2 cells. *Microbiol Immunol* 1990;34:147-51.

20. Benoit R, Dorval D, Loulergue J, et al. Post-antibiotic diarrheas: role of *Klebsiella oxytoca*. *Gastroenterol Clin Biol* 1992; 16:860-4. (In French.)

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
4. Müller H.E. 1986. Occurrence and pathogenic role of *Morganella Proteus-Providencia* group bacteria in human feces. J Clin Microbiol 23:404 – 405

5. Amin O.M. 2011. The contribution of pathogenic bacteria to GI symptoms in parasite-free patients. J Bacteriol Parasitol 2:109.

6. Mondot S, Lepage P, Seksik P, Allez M, Treton X, Bouhnik Y, Colombel J.F., Leclerc M, Pochart P, Dore J, Marteau P, GETAID. 2016. Structural robustness of the gut mucosal microbiota is associated with Crohn's disease remission after surgery. Gut 65:954 –962. <https://doi.org/10.1136/gutjnl-2015-309184>.

7. Wright E.K., Kamm M.A., Wagner J, Teo S.M., Cruz P, Hamilton A.L., Ritchie K.J., Inouye M, Kirkwood CD. 2017. Microbial factors associated with postoperative Crohn's disease recurrence. J Crohns Colitis 11:191–203.

8. Reuken P.A., Kruis W, Maaser C, Teich N, Büning J, Preiß J.C., Schmelz R, Bruns T, Fichtner-Feigl S, Stallmach A. 5 February 2018. Microbial spectrum of intra-abdominal abscesses in perforating Crohn's disease: results from a prospective German registry. J Crohns Colitis

9. Neut C, Guillemot F, Colombel J.F. 1997. Nitrate-reducing bacteria in diversion colitis: a clue to inflammation? Dig Dis Sci 42:2577–2580.

10. Neut C, Bulois P, Desreumaux P, Membree J-M, Lederman E, Gambiez L, Cortot A, Quandalle P, van Kruiningen H, Colombel J-F. 2002. Changes in the bacterial flora of the neoterminal ileum after ileocolonic resection for Crohn's disease. Am J Gastroenterol 97:939 –946

11. Hau T. 1990. Bacteria, toxins, and the peritoneum. World J Surg 14: 167–175.

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

Clostridium perfringens

1. Asha N.J., Tompkins D., Wilcox M.H. Comparative analysis of prevalence, risk factors, and molecular epidemiology of antibiotic-associated diarrhea due to *Clostridium difficile*, *Clostridium perfringens*, and *Staphylococcus aureus* // J. Clin. Microbiol. — 2006. — 44. — 2785-91.
2. Modi N., Wilcox M.H. Evidence for antibiotic induced *Clostridium perfringens* diarrhea // J. Clin. Pathol. — 2001. — 54. — 748-51.
4. Freedman J, Shrestha A, McClane B. *Clostridium perfringens* enterotoxin: action, genetics, and translational applications. *Toxins* 2016;8:73.<https://doi.org/10.3390/toxins8030073>
5. Hassan, K. A. et al. Genomic analyses of *Clostridium perfringens* isolates from five toxinotypes. *Res. Microbiol.* 166, 255–263 (2015).
6. Sim, K. et al. Dysbiosis anticipating necrotizing enterocolitis in very premature infants. *Clin. Infect. Dis.* 60, 389–397 (2015).
7. Heida, F. H. et al. A necrotizing enterocolitis-associated gut microbiota is present in the meconium: results of a prospective study. *Clin. Infect. Dis.* 62, 863–870 (2016)
8. Petit, L., Gibert, M. & Popoff, M. *Clostridium perfringens*: toxinotype and genotype. *Trends Microbiol.* 7, 104–110 (1999).

9. Kiu, R., Caim, S., Alexander, S., Pachori, P. & Hall, L. J. Probing genomic aspects of the multi-host pathogen *Clostridium perfringens* reveals significant pan-genome diversity, and a diverse array of virulence factors. *Front. Microbiol.* 8, 2485 (2017).
10. Freedman, J.C., Shrestha, A. & McClane, B.A. *Clostridium perfringens* Enterotoxin: Action, Genetics, and Translational Applications. *Toxins (Basel)* 8, <https://doi.org/10.3390/toxins8030073> (2016)
11. Freedman, J.C. et al. *Clostridium perfringens* type A-E toxin plasmids. *Res. Microbiol.* 166, 264–279 (2015).
12. Li, J. et al. Toxin plasmids of *Clostridium perfringens*. *Microbiol. Mol. Biol. Rev.* 77, 208–233 (2013).
13. Wisniewski, J. A. & Rood, J. I. The Tcp conjugation system of *Clostridium perfringens*. *Plasmid* 91, 28–36 (2017)
14. Grant, K.A. et al. The identification and characterization of *Clostridium perfringens* by real-time PCR, location of enterotoxin gene, and heat resistance. *Foodborne Pathog. Dis.* 5, 629–639 (2008).

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6. Bartlett J.G., Willey S.H., Chang T.W., Lowe B: Cephalosporin-Associated Pseudomembranous Colitis due to Clostridium difficile. *JAMA* 1979, 242:2683–5.
8. Kelly C.P., LaMont J.T. Clostridium difficile: more difficult than ever. *N Engl J Med* 2008;359:1932-1940
9. Makarewicz V.A. Clostridium difficile-associated disease in the HIV-infected patient. *J Assoc Nurses AIDS Care* 2008;19:338–41.
10. Gifford A.H, Kirkland K.B. Risk factors for Clostridium difficile-associated diarrhea on an adult hematology–oncology ward. *Eur J Clin Microbiol Infect Dis* 2006; 25:751–5.
11. Chang J.Y., Antonopoulos D.A., Kalra A, et al. Decreased diversity of the fecal microbiome in recurrent Clostridium difficile-associated diarrhea. *J Infect Dis* 2008;197:435-438

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2. Gorkiewicz G: Nosocomial and antibiotic-associated diarrhoea caused by organisms other than *Clostridium difficile*. *Int J Antimicrob Agents* 2009, 33(Suppl 1):S37–41.

3. Lin Z, Kotler D.P., Schlievert P.M., Sordillo E.M.: Staphylococcal enterocolitis: forgotten but Not gone? *Dig Dis Sci* 2010, 55:1200–7. 88.

4. Edwards L.A., O'Neill C, Furman M.A., Hicks S, Torrente F, Pérez-Machado M, Wellington E.M., Phillips A.D., Simon M.H.: Enterotoxin-producing staphylococci cause intestinal inflammation by a combination of direct epithelial cytopathy and superantigen-mediated T-cell activation. *Inflamm Bowel Dis* 2012, 18:624–40.

5. Flemming K, Ackermann G: Prevalence of enterotoxin producing *Staphylococcus aureus* in stools of patients with nosocomial diarrhea. *Infection* 2007, 35:356–8.

6. A Systematic Literature Review and Meta-analysis of Factors Associated with MRSA Colonization at Time of Hospital or ICU Admission James A. McKinnell, MD1,2,3, Loren G. Miller, MPH, MD1,2, Samantha J. Eells, MPH1,4, Eric Cui, BS5, and Susan S. Huang, MD, MPH5

7. Yang, P.-C.; Liu, T.; Wang, B.-Q.; Zhang, T.-Y.; An, Z.-Y.; Zheng, P.-Y.; Tian, D.-F. Rhinosinusitis derived Staphylococcal enterotoxin B possibly associates with pathogenesis of ulcerative colitis. *BMC Gastroenterol.* 2005, 5, doi:10.1186/1471-230X-5-28.

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2. Bisno A.L. Medical progress: group A streptococcal infections and acute rheumatic fever. *N. Engl. J. Med.* 1991;325:783–793.

3. Butler J.C., Breiman R.F., Campbell J.F. et al. Pneumococcal polysaccharide vaccine efficacy: an evaluation of current recommendations. *JAMA.* 1993;270:1826–1831.

4. CDC. Addressing emerging infectious disease threats: a prevention strategy for the United States. Atlanta: US Dept Health & Human Serv. Pub Hlth Serv. CDC 1994. CDC. Preventing the spread of vancomycin resistance — report from the Hospital Infection Control Practices Advisory Committee. *Fed Reg.* 1994; 59:25757-25763.

5. Coykendall A.L. Classification and identification of the viridans streptococci. *Clin Microbiol Rev.* 1989;2:315–328.

6. Dillon H.C. Post-streptococcal glomerulonephritis following pyoderma. *Rev Infec Dis.* 1979;1:935–943.

7. Denny F.W, Wannamaker L.W., Brink W.R. et al. Prevention of rheumatic fever: treatment of the preceding streptococcal infection. *JAMA.* 1950;143:151–153.

8. Farley M.M., Harvey R.C., Stull T. et al. A population-based assessment of invasive disease due to group B streptococcus in nonpregnant adults. *N. Engl. J. Med.* 1993;328:1807–1811.

9. Hoge C.W., Schwartz B, Talkington D.F. et al. The changing epidemiology of invasive group A streptococcal infections and the emergence of streptococcal toxic shock-like syndrome. *JAMA.* 1993;269:384–389.

10. Cartun R.W., Van Kruiningen H.J., Pedersen C.A., et al: An immunocytochemical search for infectious agents in Crohn's disease. *Mod Pathol* 1993; 6: 212–219.

11. Cucchiara, S., Iebba, V., Conte, M. P., & Schippa, S. (2009). The Microbiota in Inflammatory Bowel Disease in Different Age Groups. *Digestive Diseases*, 27(3), 252–258. doi:10.1159/000228558.

12. Thomas A.M., Manghi P, Asnicar F, Pasolli E, Armanini F, Zolfo M, et al. Metagenomic analysis of colorectal cancer datasets identifies cross-cohort microbial diagnostic signatures and a link with choline degradation. *Nat Med.* 2019;25(4):667-78.

13. Pasquereau-Kotula E, Martins M, Aymeric L, Dramsi S. Significance of *Streptococcus gallolyticus* subsp. *gallolyticus* Association With Colorectal Cancer. *Front Microbiol.* 2018;9:614.

14. Trichopoulos D, Polychronopoulou A. Epidemiology, diet and colorectal cancer. *Eur J Cancer Clin Oncol* 1986.

15. Burns C.A., McCaughey R, Lauter C.B. The association of *Streptococcus bovis* fecal carriage and colon neoplasia: possible relationship with polyps and their premalignant potential. *Am J Gastroenterol* 1985;80:42-6.

16. Boleij A, Muytjens C.M., Bukhari S.I., et al. Novel clues on the specific association of *Streptococcus gallolyticus* subsp. *gallolyticus* with colorectal cancer. *J Infect Dis* 2011.

17. Abdulmir A.S., Hafidh R.R., Mahdi L.K., et al. Investigation into the controversial association of *Streptococcus gallolyticus* with colorectal cancer and adenoma. *BMC Cancer* 2009.

18. Tjalsma H, Schöller-Guinard M, Lasonder E, et al. Profiling the humoral immune response in colon cancer patients: diagnostic antigens from *Streptococcus bovis*. *Int J Cancer* 2006.

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
1. Shiga Toxin 2-Producing *Acinetobacter haemolyticus* Associated with a Case of Bloody Diarrhea.
2. Rice L. Federal Funding for the Study of Antimicrobial Resistance in Nosocomial Pathogens: No ESKAPE. *J Infect Dis.* 2008;197(8):1079-1081. doi:10.1086/533452
3. Tang M, Poles J, Leung J et al. Inferred metagenomic comparison of mucosal and fecal microbiota from individuals undergoing routine screening colonoscopy reveals similar differences observed during active inflammation. *Gut Microbes.* 2015;6(1):48-56. doi:10.1080/19490976.2014.1000080.
4. Dominguez-Bello M, Costello E, Contreras M et al. Delivery mode shapes the acquisition and structure of the initial microbiota across multiple body habitats in newborns. *Proceedings of the National Academy of Sciences.* 2010;107(26):11971-11975. doi:10.1073/pnas.1002601107.

5. Ye L, Mueller O, Bagwell J, Bagnat M, Liddle R, Rawls J. High fat diet induces microbiota-dependent silencing of enteroendocrine cells. *Elife.* 2019;8. doi:10.7554/elife.48479.
6. Bassols J, Serino M, Carreras-Badosa G et al. Gestational diabetes is associated with changes in placental microbiota and microbiome. *Pediatr Res.* 2016;80(6):777-784. doi:10.1038/pr.2016.155.
7. Cekanaviciute E, Yoo B, Runia T et al. Gut bacteria from multiple sclerosis patients modulate human T cells and exacerbate symptoms in mouse models. *Proceedings of the National Academy of Sciences.* 2017;114(40):10713-10718. doi:10.1073/pnas.1711235114

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Prevotella spp.

1. The immune response to Prevotella bacteria in chronic inflammatory disease. Jeppe Madura Larsen.
2. Arumugam M, Raes J, Pelletier E, Le Paslier D, Yamada T, Mende DR et al Enterotypes of the human gut microbiome. Nature 2011.
3. Lozupone C.A., Li M, Campbell T.B., Flores S.C., Linderman D, Gebert M.J. et al Alterations in the gut microbiota associated with HIV – 1 infection. Cell Host Microbe 2013.
4. Ling Z, Jin C, Xie T, Cheng Y, Li L, Wu N. Alterations in the fecal microbiota of patients with HIV – 1 infection: an observational study in a Chinese population. Sci Rep 2016.
5. Yang L, Poles M.A., Fisch G.S., Ma Y, Nossa C, Phelan J.A. et al HIV – induced immunosuppression is associated with colonization of the proximal gut by environmental bacteria. AIDS 2016.
6. Dillon S.M., Lee E.J., Kotter C.V., Austin G.L., Dong Z, Hecht D.K. et al An altered intestinal mucosal microbiome in HIV – 1 infection is associated with mucosal and systemic immune activation and endotoxemia. Mucosal Immunol 2014.
7. Dillon S.M., Lee E.J., Kotter C.V., Austin G.L., Gianella S, Siewe B et al Gut dendritic cell activation links an altered colonic microbiome to mucosal and systemic T-cell activation in untreated HIV – 1 infection. Mucosal Immunol 2016.
8. Pedersen H.K., Gudmundsdottir V, Nielsen H.B., Hyotylainen T, Nielsen T, Jensen BAH et al Human gut microbes impact host serum metabolome and insulin sensitivity. Nature 2016.

9. Moreno – Indias I, Sánchez – Alcoholado L, García – Fuentes E, Cardona F, Queipo – Ortuño M.I., Tinahones F.J. Insulin resistance is associated with specific gut microbiota in appendix samples from morbidly obese patients. Am J Transl Res 2016.
10. Zhu L, Baker S.S., Gill C, Liu W, Alkhouri R, Baker R.D. et al Characterization of gut microbiomes in nonalcoholic steatohepatitis (NASH) patients: a connection between endogenous alcohol and NASH. Hepatology 2013.
11. Hu H – J, Park S – G, Jang H.B., Choi M – G, Park K – H, Kang J.H. et al Obesity alters the microbial community profile in Korean adolescents. PLoS One 2015.
12. Li J, Zhao F, Wang Y, Chen J, Tao J, Tian G et al Gut microbiota dysbiosis contributes to the development of hypertension. Microbiome 2017.
13. Michail S, Lin M, Frey M.R., Fanter R, Paliy O, Hilbush B et al Altered gut microbial energy and metabolism in children with non-alcoholic fatty liver disease. FEMS Microbiol Ecol 2015.
14. Forslund K, Hildebrand F, Nielsen T, Falony G, Le Chatelier E, Sunagawa S et al Disentangling type 2 diabetes and metformin treatment signatures in the human gut microbiota. Nature 2015.

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
4. Avivi-Green, C., Polak-Charcon, S., Madar, Z. & Schwartz, B. (2000;). Apoptosis cascade proteins are regulated in vivo by high intracolonic butyrate concentration: correlation with colon cancer inhibition. *Oncol Res* 12, 83–95.

5. Pryde, S. E., Duncan, S. H., Hold, G. L., Stewart, C. S. & Flint, H. J. (2002;). The microbiology of butyrate formation in the human colon. *FEMS Microbiol Lett* 217, 133–139.

6. Knoll R.L., Forslund K, Kultima J.R., Meyer C.U., Kullmer U, Sunagawa S, et al. Gut microbiota differs between children with inflammatory bowel disease and healthy siblings in taxonomic and functional composition: a metagenomic analysis. *Am J Physiol Gastrointest Liver Physiol* (2017) 312:G327–39. doi:10.1152/ajpgi.00293.2016.

7. Vermeiren J, Van den Abbeele P, Laukens D, Vigsnaes L.K., De Vos M, Boon N, et al. Decreased colonization of fecal *Clostridium coccoides*/*Eubacterium rectale* species from ulcerative colitis patients in an in vitro dynamic gut model with mucin environment. *FEMS Microbiol Ecol* (2012)

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2. Ruminococcus gnavus, a member of the human gut microbiome associated with Crohn's disease, produces an inflammatory polysaccharide. Henke M.T., Kenny D.J., Cassilly C.D., Vlamakis H, Xavier R.J., Clardy J.

3. Intestinal Dysbiosis Featuring Abundance of Ruminococcus gnavus Associates With Allergic Diseases in Infants. Chua H.H., Chou H.C., Tung Y.L., Chiang B.L., Liao C.C., Liu H.H., Ni Y.H.

4. j Nagao-Kitamoto H, Kamada N (February 2017). «Host-microbial Cross-talk in Inflammatory Bowel Disease». Immune Network. 17 (1): 1–12. doi:10.4110/in.2017.17.1.

5. Hansen et al., 2013.

6. Titécat et al. 2014.

7. Vereecke and Elewaut, 2017.

8. Eckburg et al., 2005; Lay et al., 2005.

9. Kang et al., 2010.

10. Qin et al. al., 2010

1. Barcenilla, S.E. Pryde, J.C. Martin, S.H. Duncan, C.S. Stewart, C. Henderson, H.J. Flint, Phylogenetic relationships of butyrate-producing bacteria from the human gut, Appl. Environ. Microbiol. 66 (2000) 1654e1661.

2. Reduced Abundance of Butyrate-Producing Bacteria Species in the Fecal Microbial Community in Crohn's Disease. Takahashi K, Nishida A, Fujimoto T, Fujii M, Shioya M, Imaeda H, Inatomi O, Bamba S, Sugimoto M, Andoh A.

3. Gut Microbiome Function Predicts Response to Anti-integrin Biologic Therapy in Inflammatory Bowel Diseases. Ananthakrishnan A.N., Luo C, Yajnik V, Khalili H, Garber J.J., Stevens B.W., Cleland T, Xavier R.J.

4. Microbiota in anorexia nervosa: The triangle between bacterial species, metabolites and psychological tests. Borgo F, Riva A, Benetti A, Casiraghi M.C., Bertelli S, Garbossa S, Anselmetti S, Scarone S, Pontiroli A.E., Morace G, Borghi E.

5. Gut Dysbiosis in Patients with Anorexia Nervosa. Morita C, Tsuji H, Hata T, Gondo M, Takakura S, Kawai K, Yoshihara K, Ogata K, Nomoto K, Miyazaki K, Sudo N PLoS One. 2015;

6. Weight gain in anorexia nervosa does not ameliorate the faecal microbiota, branched chain fatty acid profiles, and gastrointestinal complaints. Mack I, Cuntz U, Grämer C, Niedermaier S, Pohl C, Schwiertz A, Zimmermann K, Zipfel S, Enck P, Penders J Sci Rep. 2016 May 27

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1. Ansel Hsiao, Jun Zhu. Pathogenicity and virulence regulation of *Vibrio cholerae* at the interface of host-gut microbiome interactions. *Virulence*. 2020 Dec;11(1):1582-1599. doi: 10.1080/21505594.2020.1845039.

2. Jenq R.R., Taur Y, Devlin S.M., Ponce D.M., Goldberg J.D., Ahr K.F., Littmann E.R., Ling L, Gobourne A.C., Miller L.C., Docampo M.D., Peled J.U., Arpaia N, Cross J.R., Peets T.K., Lumish M.A., Shono Y, Dudakov J.A., Poeck H, Hanash A.M., Barker J.N., Perales M-A, Giral S.A., Pamer E.G., van den Brink MRM. Intestinal *Blautia* Is Associated with Reduced Death from Graft-versus-Host Disease. *Biol Blood Marrow Transplant*. 2015 Aug; 21(8):1373-83. doi: 10.1016/j.bbmt.2015.04.016.

3. Scheithauer TPM, Rampanelli E, Nieuwdorp M, Vallance B.A., Verchere C.B., van Raalte D.H.,

Herrema H. Gut Microbiota as a Trigger for Metabolic Inflammation in Obesity and Type 2 Diabetes. *Front Immunol*. 2020 Oct 16;11:571731. doi: 10.3389/fimmu.2020.571731.

4. Pedersen H.K., Gudmundsdottir V, Nielsen H.B., Hyotylainen T, Nielsen T, Jensen B.A., et al. Human gut microbes impact host serum metabolome and insulin sensitivity. *Nature* (2016) 535(7612):376–81. doi: 10.1038/nature18646.

5. Yuan X, Chen R, Zhang Y, Lin X, Yang X, McCormick K.L. Gut Microbiota of Chinese Obese Children and Adolescents With and Without Insulin Resistance. *Front Endocrinol (Lausanne)*. 2021 Mar 19;12:636272. doi: 10.3389/fendo.2021.636272.

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
2. Chaudhary P.P., Sirohi S.K., Saxena J (2012) Diversity analysis of methanogens in rumen of *Bubalus bubalis* by 16S ribotyping and sequence analysis. *Gene* 493:13–17. <https://doi.org/10.1016/j.gene.2011.11.041>

3. Lecours P.B., Marsolais D, Cormier Y, Berberi M, Haché C, Bourdages R, Duchaine C (2014) Increased prevalence of *Methanospiraeta* in inflammatory bowel diseases. *PLoS One* 9:e87734. <https://doi.org/10.1371/journal.pone.0087734>

4. Pimentel M, Gunsalus R.P., Rao S.S., Zhang H (2012) Methanogens in human health and disease. *Am J Gastroenterol Suppl* 1:28–33. <https://doi.org/10.1038/ajgsup.2012.6>

5. Pimentel M, Mayer A.G., Park S, Chow E.J., Hasan A, Kong Y (2003) Methane production during lactulose breath test is associated with gastrointestinal disease presentation. *Dig Dis Sci* 48:86–92. <https://doi.org/10.1023/A:1021738515885>

6. Lais-Lecours P, Marsolais D, Cormier Y, Berberi M, Hache C, et al. (2014) Increased Prevalence of *Methanospiraeta* in Inflammatory Bowel Diseases. *PLOS ONE* 9: e87734.

7. Mathur R, Goyal D, Kim G, Barlow G.M., Chua K.S., Pimentel M. Methane-producing human subjects have higher serum glucose levels during oral glucose challenge than non-methane producers: a pilot study of the effects of enteric methanogens on glycemic regulation. *Research Journal of Endocrinology and Metabolism*. 2014;2:1-5.

8. Cesario V, Di Rienzo T.A., Campanale M, et al. Methane intestinal production and poor metabolic control in type I diabetes complicated by autonomic neuropathy. *Minerva Endocrinol*. 2014;39:201-207.

Armougom F., Henry M, Vialettes B, Raccach D, Raoult D. Monitoring bacterial community of human gut microbiota reveals an increase in *Lactobacillus* in obese patients and Methanogens in anorexic patients. *PLoS One*. 2009 Sep 23;4(9):e7125. doi: 10.1371/journal.pone.0007125.

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2. Purcell R.V., Visnovska M, Biggs P.J., Schmeier S, Frizelle F.A. 2017. Distinct gut microbiome patterns associate with consensus molecular subtypes of colorectal cancer. *Sci Rep* 7:11590. <https://doi.org/10.1038/s41598-017-11237-6>.

3. Feng Q, Liang S.S., Jia H.J., Stadlmayr A, Tang L.Q., Lan Z, Zhang D.Y., Xia H.H., Xu X.Y., Jie Z.Y., et al. Gut microbiome development along the colorectal adenoma-carcinoma sequence. *Nat Commun*. 2015;6

4. Gao Z, Guo B, Gao R, Zhu Q, Qin H. Microbiota dysbiosis is associated with colorectal cancer. *Front Microbiol*. 2015;6:20.

5. Yu J, Feng Q, Wong S.H., Zhang D, Yi Liang Q, Qin Y, Tang L, Zhao H, Stenvang J, Li Y. Metagenomic analysis of faecal microbiome as a tool towards targeted non-invasive biomarkers for colorectal cancer. *Gut*. 2015; gutjnl-2015-309800

6. Distinct gut microbiome patterns associate with consensus molecular subtypes of colorectal cancer. Purcell R.V., Visnovska M, Biggs P.J., Schmeier S, Frizelle F.A.

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
Fusobacterium nucleatum

1. Gevers D et al. The treatment-naive microbiome in new-onset Crohn's disease. *Cell Host Microbe* 15, 382–392 (2014). [PubMed: 24629344]
2. Pascal V et al. A microbial signature for Crohn's disease. *Gut* 66, 813–822 (2017). [PubMed: 28179361]
3. Strauss J et al. Invasive potential of gut mucosa-derived *Fusobacterium nucleatum* positively correlates with IBD status of the host. *Inflamm. Bowel Dis* 17, 1971–1978 (2011). [PubMed: 21830275]
4. McCoy A.N. et al. *Fusobacterium* is associated with colorectal adenomas. *PLOS ONE* 8, e53653 (2013). [PubMed: 23335968]
5. Kostic A.D. et al. Genomic analysis identifies association of *Fusobacterium* with colorectal carcinoma. *Genome Res* 22, 292–298 (2012). [PubMed: 22009990]
6. Castellarin M et al. *Fusobacterium nucleatum* infection is prevalent in human colorectal carcinoma. *Genome Res* 22, 299–306 (2012). [PubMed: 22009989]
7. Feng Q et al. Gut microbiome development along the colorectal adenoma-carcinoma sequence. *Nat. Commun* 6, 6528 (2015). [PubMed: 25758642]
8. Yu J et al. Invasive *Fusobacterium nucleatum* may play a role in the carcinogenesis of proximal colon cancer through the serrated neoplasia pathway. *Int. J. Cancer* 139, 1318–1326 (2016). [PubMed: 27130618]

9. Guo S et al. A simple and novel fecal biomarker for colorectal cancer: ratio of *Fusobacterium nucleatum* to probiotics populations, based on their antagonistic effect. *Clin. Chem* 64, 1327–1337 (2018). [PubMed: 29914865]
10. Mima K et al. *Fusobacterium nucleatum* in colorectal carcinoma tissue and patient prognosis. *Gut* 65, 1973–1980 (2015). [PubMed: 26311717] This molecular epidemiology study highlights the association of tumoural fusobacterial levels with worse prognosis in colorectal cancer.
11. Wang H-F et al. Evaluation of antibody level against *Fusobacterium nucleatum* in the serological diagnosis of colorectal cancer. *Sci. Rep* 6, 33440 (2016). [PubMed: 27678333]
12. Kostic A.D., Chun E, Robertson L, et al. *Fusobacterium nucleatum* potentiates intestinal tumorigenesis and modulates the tumor-immune microenvironment. *Cell Host Microbe* 2013; 14:207–15.
13. Rubinstein M.R., Wang X, Liu W, et al. *Fusobacterium nucleatum* promotes colorectal carcinogenesis by modulating E-cadherin/beta-catenin signaling via its FadA adhesin. *Cell Host Microb* 2013;14:195–206.
14. Abed J, Emgård J.E., Zamir G, et al. Fap2 mediates *Fusobacterium nucleatum* colorectal adenocarcinoma enrichment by binding to tumor-expressed Gal-GalNAc. *Cell Host Microbe* 2016;20:215–25

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